


# Other rhythm


Provisional on-device interpretation

| RECORDED             | DURATION   | QUALITY | CONFIDENCE |
|----------------------|------------|---------|------------|
| 15 Jul 2025 at 18:27 | 30 seconds | 94%     | 100%       |


MODEL VIEW



30-SECOND RHYTHM STRIP



25 mm/s / 10 mm/mV / 256 Hz



HOW TO READ

Each row is five seconds at a true 25 mm/s time base. The same vertical gain is used throughout. Model influence shows which masked interval most changed the winning score.

READ WITH CARE

This research-only software result is not a diagnosis. Seek urgent care for chest pain, fainting, severe breathlessness, or other concerning symptoms. A clinician should interpret this single-lead recording in context.