

Atrial fibrillation

Provisional on-device interpretation


RECORDED	DURATION	QUALITY	CONFIDENCE
15 Jul 2025, 18:24	30 seconds	94%	100%

MODEL VIEW



30-SECOND RHYTHM STRIP

25 mm/s / 10 mm/mV / 256 Hz



HOW TO READ

Each row is five seconds at a true 25 mm/s time base. The same vertical gain is used throughout. Model influence shows which masked interval most changed the winning score.

READ WITH CARE

This research-only software result is not a diagnosis. Seek urgent care for chest pain, fainting, severe breathlessness, or other concerning symptoms. A clinician should interpret this single-lead recording in context.